



NORSE RESCUE® NASAL AIRWAY NPA 28F 9.3 MM LATEX FREE, STERILE



Item No.:	AW-00057-1
Brand:	NORSE RESCUE®
Manufacturer:	EMSRUN GROUP
Certificates:	CE, MDR, FDA
Classification:	Class Is
Size:	28F / 9.3 mm
Color:	Green
Country of origin:	CN
NSN:	6515-22-640-9390
Tariff:	90192090
Sales Units:	1 pc, 20 pcs, 400 pcs
Weight:	10 g (per unit)
Weight (Gross):	11
Dimensions:	17 × 6 × 2.5 cm
Shelf life:	4 Yrs from prod. date

The NORSE RESCUE® Nasal Airway (NPA) size 28F is an essential emergency tool, designed for patient comfort and optimal airflow. Made from soft PVC, it minimizes discomfort and is engineered to maintain shape, preventing collapse and ensuring efficient oxygenation and ventilation in critical situations.

AP Services A/S

Niels P. Thomsens Vej 3, 7500 Holstebro, Denmark
+45 96 90 70 90 | ap-services.dk

