

NORSE RESCUE® GUEDEL, SIZE 4, RED, 100 MM



Item No.:	AW-00120
Brand:	NORSE RESCUE®
Manufacturer:	EMSRUN GROUP
Certificates:	CE, MDR, FDA
Classification:	Class I
Size:	Size 4 (ISO10)
Color:	Red
Country of origin:	CN
NSN:	6515-22-640-9391
Tariff:	90183900
Sales Units:	1 pc, 250 pcs
Weight:	17 g (per unit)
Weight (Gross):	20
Dimensions:	11 × 4.5 × 4 cm
Shelf life:	N/A

The NORSE RESCUE® Oropharyngeal Airway (OPA) is a critical tool for maintaining unobstructed airflow in unconscious patients, preventing tongue-related airway blockage. It features a central lumen for suction catheter insertion, ensuring effective secretion clearance and optimal patient ventilation and oxygenation.

