

NORSE RESCUE® RESCUE BREATHER BLACK HARD CASE INCL. GLOVES



Item No.:	BR-0011
Brand:	NORSE RESCUE®
Manufacturer:	EMSRUN GROUP
Certificates:	CE, MDR, FDA
Classification:	Class I
Size:	Adult
Color:	Transparent, Black
Country of origin:	CN
NSN:	4240-22-632-9655
Tariff:	90192090
Sales Units:	1 pc, 60 pcs
Weight:	122 g (per unit)
Weight (Gross):	136
Dimensions:	15 × 11.5 × 5 cm
Shelf life:	5 Yrs from prod. date

The NORSE RESCUE® Rescue Breather is an emergency tool designed for effective mouth-to-mouth ventilation, featuring an oxygen inlet for enhanced support. Housed in a tactical black clamshell case for rapid accessibility, this high-quality mask is accompanied by gloves to ensure hygiene and protection during life-saving procedures.

