

NORSE RESCUE® CHEST SEAL, VENTED 1 PC - FOUR CHANNEL, W/ GAUZE



Item No.:	BR-00187
Brand:	NORSE RESCUE®
Manufacturer:	EMSRUN GROUP
Certificates:	CE, MDR, FDA
Classification:	Class I
Size:	15 x 15 cm
Color:	Transparent
Country of origin:	CN
NSN:	6510-22-636-8498
Tariff:	30051000
Sales Units:	1 pc, 20 pcs, 200 pcs
Weight:	37 g (per unit)
Weight (Gross):	37
Dimensions:	14 x 10 x 1 cm
Shelf life:	5 Yrs from prod. date

The NORSE RESCUE® Chest Seal Vented is a critical medical dressing for open chest wounds, designed to prevent tension pneumothorax by maintaining a secure seal under extreme conditions. Its four-vented design allows air to escape while preventing air inflow, ensuring optimal patient care in traumatic situations.

AP Services A/S

Niels P. Thomsens Vej 3, 7500 Holstebro, Denmark
+45 96 90 70 90 | ap-services.dk

