

NORSE RESCUE® CHEST SEAL, VENTED 2PCS - FOUR CHANNEL, W/GAUZE



Item No.:	BR-00187-TWIN
Brand:	NORSE RESCUE®
Manufacturer:	WUXI EMSRUN
Certificates:	CE, MDR
Classification:	Class I
Size:	15 x 15 cm
Color:	Transparent
Country of origin:	CN
Tariff:	30051000
Sales Units:	1 pc
Weight:	90 g (per unit)
Weight (Gross):	97
Dimensions:	15 x 11 x 2 cm
Shelf life:	5 Yrs from prod. date

The NORSE RESCUE® Chest Seal Four Channel Vented - Twin Pack is an advanced occlusive dressing designed to treat open chest wounds effectively. Equipped with a medical-grade hydrogel adhesive, it ensures secure adhesion even on skin exposed to blood, sweat, hair, sand, or water. Each pack includes gauze for additional wound care.

AP Services A/S

Niels P. Thomsens Vej 3, 7500 Holstebro, Denmark
+45 96 90 70 90 | ap-services.dk

