

NORSE RESCUE® CHEST SEAL

1 PC W/GAUZE, 15 X 15 CM



Item No.:	BR-00188
Brand:	NORSE RESCUE®
Manufacturer:	EMSRUN GROUP
Certificates:	CE, MDR, FDA
Classification:	Class I
Size:	15 x 15 cm
Color:	Transparent
Country of origin:	CN
NSN:	6510-22-636-8506
Tariff:	30051000
Sales Units:	1 pc, 20 pcs, 200 pcs
Weight:	42 g (per unit)
Weight (Gross):	42
Dimensions:	14 x 10 x 1 cm
Shelf life:	5 Yrs from prod. date

The NORSE RESCUE® Chest Seal Non-Vented is a specially designed advanced dressing for the treatment of open chest wounds, aiming to prevent tension pneumothorax. Crafted with a medical-grade hydrogel adhesive, it ensures a robust seal across various conditions, providing reliable wound management and enhancing patient safety in traumatic situations.

