

## NORSE RESCUE® EMERGENCY CHEST DRAIN KIT



<b>Item No.:</b>	BR-100060
<b>Brand:</b>	NORSE RESCUE®
<b>Procedure Pack</b>	
<b>Producer:</b>	AP Services A/S
<b>Country of origin:</b>	DK
<b>NSN:</b>	6515-22-639-2163
<b>Tariff:</b>	90189084
<b>Sales Units:</b>	1 pc
<b>Weight:</b>	214 g (per unit)
<b>Weight (Gross):</b>	214
<b>Dimensions:</b>	30 x 18 x 3 cm
<b>Shelf life:</b>	N/A

The NORSE RESCUE® Emergency Chest Drainage Kit is an all-encompassing solution for immediate and effective chest drainage procedures. This kit contains essential, sterile equipment including a CH32 drainage set with a 600 ml collection bag, a safety scalpel, suture materials, durable nitrile gloves, a sharp Pean clamp, disinfection wipes, and a specialized chest seal to prevent air from entering the pleural space. Designed for quick response and safety, this kit is indispensable for emergency medical situations.

