

NORSE RESCUE® OXYGEN MASK, PEDIATRIC PEDIATRIC, 200 CM HOSE



| | |
|---------------------------|-----------------------|
| Item No.: | BR-100077 |
| Brand: | NORSE RESCUE® |
| Manufacturer: | Ningbo MFLAB Medical |
| Certificates: | CE, MDD |
| Classification: | Class IIa |
| Size: | M, Pediatric |
| Color: | Transparent |
| Country of origin: | CN |
| NSN: | 4240-22-640-9426 |
| Tariff: | 90192090 |
| Sales Units: | 1 pc, 100 pcs |
| Weight: | 54 g (per unit) |
| Weight (Gross): | 56 |
| Dimensions: | 16 × 19 × 3 cm |
| Shelf life: | 5 Yrs from prod. date |

The NORSE RESCUE® Oxygen Mask, Medium (Pediatric), is designed for children requiring medium concentration oxygen therapy. Ideal for situations such as awakening, minor trauma, and patient transport, this mask delivers oxygen concentrations between 40-60%, depending on breathing patterns, with a hose length of 200 cm for mobility and comfort.

