

NORSE RESCUE® OXYGEN MASK W/ RESERVOIR PE CHILD, THREE IN ONE EMERGENCY MASK



Item No.:	BR-100078
Brand:	NORSE RESCUE®
Manufacturer:	Ningbo MFLAB Medical
Certificates:	CE, MDD
Classification:	Class IIa
Size:	Medium, Pediatric
Color:	Transparent
Country of origin:	CN
NSN:	4240-22-640-9458
Tariff:	90192090
Sales Units:	1 pc, 100 pcs
Weight:	63 g (per unit)
Weight (Gross):	65
Dimensions:	16 × 19 × 3 cm
Shelf life:	5 Yrs from prod. date

The NORSE RESCUE® Oxygen Mask with Reservoir, size Medium, is a pediatric three-in-one emergency mask designed for flexibility in oxygen delivery. It can be adjusted from high to medium concentration levels, equipped with a side valve, reservoir bag, and a standard connector, complemented by a 200 cm hose for extensive mobility and ease of use.

