

NORSE RESCUE® OXYGEN MASK W/RESERVOIR ADU ADULT, THREE IN ONE EMERGENCY MASK



Item No.:	BR-100079
Brand:	NORSE RESCUE®
Manufacturer:	Ningbo MFLAB Medical
Certificates:	CE, MDD
Classification:	Class IIa
Size:	L3, Adult
Color:	Transparent
Country of origin:	CN
NSN:	4240-22-640-9459
Tariff:	90192090
Sales Units:	1 pc, 100 pcs
Weight:	75 g (per unit)
Weight (Gross):	77
Dimensions:	18 × 21 × 4 cm
Shelf life:	5 Yrs from prod. date

The NORSE RESCUE® Oxygen Mask w/Reservoir L3 for adults is a versatile three-in-one emergency mask designed to deliver oxygen treatment at approximately 10 liters per minute. Featuring a soft construction with a metal nose clip for a secure fit and two side vents for optimal oxygen concentration, this mask ensures effective treatment with a 200 cm hose for extended reach.

