

NORSE RESCUE® NEBULIZER MASK, ADULT ADULT, 200 CM HOSE



Item No.:	BR-100081
Brand:	NORSE RESCUE®
Manufacturer:	Ningbo MFLAB Medical
Certificates:	CE, MDD
Classification:	Class IIa
Size:	L3, Adult
Color:	Transparent
Country of origin:	CN
Tariff:	90192090
Sales Units:	1 pc, 100 pcs
Weight:	78 g (per unit)
Weight (Gross):	80
Dimensions:	18 × 24 × 4 cm
Shelf life:	5 Yrs from prod. date

The NORSE RESCUE® Nebulizer Mask L3 for adults is designed to provide effective nebulization treatment, suitable for oxygen supplementation and medication administration. It ensures efficient medication absorption through nebulization, equipped with a 200 cm hose for enhanced mobility and comfort during use.

