

NORSE RESCUE® DECOPIN 10GA DECOMPRESSION NEEDLE



Item No.:	BR-200010
Brand:	NORSE RESCUE®
Manufacturer:	EMSRUN GROUP
Certificates:	CE, MDR
Classification:	Class IIa
Size:	10ga 3.75"
Color:	Grey, Red
Country of origin:	CN
NSN:	4240-22-640-9479
Tariff:	90183900
Sales Units:	1 pc, 10 pcs, 200 pcs
Weight:	28 g (per unit)
Weight (Gross):	31
Dimensions:	20 × 2.5 × 2 cm
Shelf life:	5 Yrs from prod. date

The NORSE RESCUE® DECOPIN is a sterile, single-use decompression needle designed for trained personnel to treat tension pneumothorax in adult trauma patients. Featuring a 10-gauge, 3.75-inch cannula, it enables rapid pressure relief in life-threatening situations and is sterilized by irradiation for safe, effective use in the field.

