

NORSE RESCUE® DECOPIN 14GA DECOMPRESSION NEEDLE



Item No.:	BR-200011
Brand:	NORSE RESCUE®
Manufacturer:	EMSRUN GROUP
Certificates:	CE, MDR
Classification:	Class IIa
Size:	14ga 3.25"
Color:	Grey, Red
Country of origin:	CN
NSN:	4240-22-640-9480
Tariff:	90183900
Sales Units:	1 pc, 10 pcs, 200 pcs
Weight:	20 g (per unit)
Weight (Gross):	22
Dimensions:	15 × 2.2 × 1.7 cm
Shelf life:	5 Yrs from prod. date

The NORSE RESCUE® DECOPIN is a sterile, single-use decompression needle designed for trained personnel to alleviate tension pneumothorax in adult patients experiencing progressive respiratory distress due to thoracic trauma. Featuring a 14-gauge, 3.25-inch cannula, it enables rapid pressure relief in life-threatening situations and is sterilized by irradiation for safe, effective use in the field.

