

NORSE RESCUE® TRAUMA BANDAGE 8" ABDOMINAL



Item No.:	CL-00458
Brand:	NORSE RESCUE®
Manufacturer:	EMSRUN GROUP
Certificates:	CE, MDR, FDA
Classification:	Class Is
Size:	8" / 20 cm, Abdominal
Color:	Tan
Country of origin:	CN
NSN:	6510-22-640-9339
Tariff:	30059050
Sales Units:	1 pc, 70 pcs
Weight:	156 g (per unit)
Weight (Gross):	168
Dimensions:	21 x 10 x 5 cm
Shelf life:	5 Yrs from prod. date

The NORSE RESCUE® Trauma Bandage 8" Abdominal is engineered for severe injuries such as abdominal eviscerations, traumatic amputations, and large wounds. Featuring an 8" elastic wrap and a large 12" by 12" sterile pad, it offers superior coverage and support. The integrated quick grip roll control ensures efficient application, complemented by an occlusive layer for enhanced abdominal wound care.

