



AEROcase® – RBVM1C | AEROcase® – SEG XL

Dear Distribution Partners,

We are glad to inform you that we are changing our **AEROcase® – RBVM1C Emergency Bag** of the 1R series from polyester to the extra robust tarpaulin material with immediate effect. Also, the plastic insert will be omitted and replaced by a tarpaulin division! We would also like to introduce the **AEROcase® – SEG Emergency Backpack**, which will soon expand our EMS series!

AEROcase® – RBVM1C of the **1R** Series - Soon available in Tarpaulin



HT03-RBVM1C-R

(Signal Red - RAL 3001)

Dimensions: approx. 47 x 28 x 25 cm

Weight: approx. 2.48 kg (incl. modular bags)

Volume: approx. 33 l



Description:

Extra robust emergency bag in "Medium" made of AEROTex® – PLAN (Tarpaulin), as well as integrated interior division made of tarpaulin material, with 3 removable partitions for an individual compartment separation. Additionally with integrated holder for oxygen cylinder 0.8 | 1 | 2 liters. Delivery including 4 modular bags.

Special features at a glance

- 4-colour coded modular bags with carrying straps for the individual organization - visible, portable and stable
- Inner compartment in the inside modular bag cover / closure flap
- Large transparent compartment on the front for individual labelling
- Equipped with carrying handles, a combined removable shoulder / back strap and 5 plastic feet
- Bidirectional zippers
- Available in Signal Red